

Package ‘lievens’

May 3, 2024

Title Real-Time PCR Data Sets by Lievens et al. (2012)

Version 0.0.1

Description Real-time quantitative polymerase chain reaction (qPCR) data sets by Lievens et al. (2012) <[doi:10.1093/nar/gkr775](https://doi.org/10.1093/nar/gkr775)>. Provides one single tabular tidy data set in long format, encompassing three dilution series, targeted against the soybean Lectin endogene. Each dilution series was assayed in one of the following PCR-efficiency-modifying conditions: no PCR inhibition, inhibition by isopropanol and inhibition by tannic acid. The inhibitors were co-diluted along with the dilution series. The co-dilution series consists of a five-point, five-fold serial dilution. For each concentration there are 18 replicates. Each amplification curve is 60 cycles long. Original raw data file is available at the Supplementary Data section at Nucleic Acids Research Online <[doi:10.1093/nar/gkr775](https://doi.org/10.1093/nar/gkr775)>.

License CC BY 4.0

Encoding UTF-8

RoxxygenNote 7.3.1

Imports tibble

Depends R (>= 2.10)

LazyData true

URL <https://rmagno.eu/lievens/>, <https://github.com/ramiromagno/lievens>

BugReports <https://github.com/ramiromagno/lievens/issues>

NeedsCompilation no

Author Ramiro Magno [aut, cre] (<<https://orcid.org/0000-0001-5226-3441>>),
Pattern Institute [cph, fnd]

Maintainer Ramiro Magno <rmagno@pattern.institute>

Repository CRAN

Date/Publication 2024-05-03 13:10:02 UTC

R topics documented:

lievens	2
---------	-------	---

lievens

*qPCR data sets by Lievens et al. (2012)***Description**

One single tabular tidy data set in long format, encompassing three data sets of five-point, five-fold dilution series: (i) without any inhibitor, (ii) with isopropanol inhibition and (iii) with tannic acid inhibition. The target amplicon consisted of a sequence within the soybean Lectin endogene. Please read the Methods section of Lievens et al. (2012) for more experimental details.

Each data set consists of a five-point, five-fold dilution series spanning an amplicon copy number range from 100,000 down to 160. Each concentration is replicated 18 times. Each reaction has been amplified through 60 cycles.

Dilution series:

```
dplyr::filter(lievens, inhibitor == "none")
#> # A tibble: 5,400 x 13
#>   plate well  target dye  inhibitor inhibitor_conc sample sample_type
#>   <fct> <fct> <fct> <fct> <fct>           <dbl> <fct> <fct>
#> 1 soy    <NA>  Le1   SYBR  none            0     S1    std
#> 2 soy    <NA>  Le1   SYBR  none            0     S1    std
#> 3 soy    <NA>  Le1   SYBR  none            0     S1    std
#> 4 soy    <NA>  Le1   SYBR  none            0     S1    std
#> 5 soy    <NA>  Le1   SYBR  none            0     S1    std
#> 6 soy    <NA>  Le1   SYBR  none            0     S1    std
#> 7 soy    <NA>  Le1   SYBR  none            0     S1    std
#> 8 soy    <NA>  Le1   SYBR  none            0     S1    std
#> 9 soy    <NA>  Le1   SYBR  none            0     S1    std
#> 10 soy   <NA>  Le1   SYBR  none            0     S1    std
#> # i 5,390 more rows
#> # i 5 more variables: replicate <fct>, copies <int>, dilution <int>,
#> #   cycle <int>, fluor <dbl>
```

Isopropanol inhibition:

A series of reactions subjected to inhibition by isopropanol with concentrations: 2.5, 0.5, 0.1, 0.02, and 0.004 % (v/v). Because samples have been co-diluted, the initial copy numbers of the target amplicon also follow the same five-fold progression in tandem: 100,000, 20,000, 4,000, 800 and 160 copies.

```
dplyr::filter(lievens, inhibitor == "isopropanol")
#> # A tibble: 5,400 x 13
#>   plate      well  target dye  inhibitor inhibitor_conc sample sample_type
#>   <fct>      <fct> <fct> <fct> <fct>           <dbl> <fct> <fct>
#> 1 soy+isopropan~ <NA>  Le1   SYBR  isopropa~        2.5  S1    std
#> 2 soy+isopropan~ <NA>  Le1   SYBR  isopropa~        2.5  S1    std
#> 3 soy+isopropan~ <NA>  Le1   SYBR  isopropa~        2.5  S1    std
```

```
#> 4 soy+isopropan~ <NA> Le1 SYBR isopropa~ 2.5 S1 std
#> 5 soy+isopropan~ <NA> Le1 SYBR isopropa~ 2.5 S1 std
#> 6 soy+isopropan~ <NA> Le1 SYBR isopropa~ 2.5 S1 std
#> 7 soy+isopropan~ <NA> Le1 SYBR isopropa~ 2.5 S1 std
#> 8 soy+isopropan~ <NA> Le1 SYBR isopropa~ 2.5 S1 std
#> 9 soy+isopropan~ <NA> Le1 SYBR isopropa~ 2.5 S1 std
#> 10 soy+isopropan~ <NA> Le1 SYBR isopropa~ 2.5 S1 std
#> # i 5,390 more rows
#> # i 5 more variables: replicate <fct>, copies <int>, dilution <int>,
#> #   cycle <int>, fluor <dbl>
```

Tannic acid inhibition:

A series of reactions subjected to inhibition by tannic acid with concentrations: 0.2, 0.04, 0.008, 0.0016 and 0.0032 μ L/mL. Because samples have been co-diluted, the initial copy numbers of the target amplicon also follow the same five-fold progression in tandem: 100,000, 20,000, 4,000, 800 and 160.

```
dplyr::filter(lievens, inhibitor == "tannic acid")
#> # A tibble: 5,400 x 13
#>   plate      well target dye  inhibitor inhibitor_conc sample sample_type
#>   <fct>     <fct> <fct> <fct> <fct> <dbl> <fct> <fct>
#> 1 soy+tannic ac~ <NA> Le1  SYBR  tannic a~    0.2 S1 std
#> 2 soy+tannic ac~ <NA> Le1  SYBR  tannic a~    0.2 S1 std
#> 3 soy+tannic ac~ <NA> Le1  SYBR  tannic a~    0.2 S1 std
#> 4 soy+tannic ac~ <NA> Le1  SYBR  tannic a~    0.2 S1 std
#> 5 soy+tannic ac~ <NA> Le1  SYBR  tannic a~    0.2 S1 std
#> 6 soy+tannic ac~ <NA> Le1  SYBR  tannic a~    0.2 S1 std
#> 7 soy+tannic ac~ <NA> Le1  SYBR  tannic a~    0.2 S1 std
#> 8 soy+tannic ac~ <NA> Le1  SYBR  tannic a~    0.2 S1 std
#> 9 soy+tannic ac~ <NA> Le1  SYBR  tannic a~    0.2 S1 std
#> 10 soy+tannic ac~ <NA> Le1  SYBR  tannic a~   0.2 S1 std
#> # i 5,390 more rows
#> # i 5 more variables: replicate <fct>, copies <int>, dilution <int>,
#> #   cycle <int>, fluor <dbl>
```

Format

A [tibble](#) providing amplification curve data in long format. Each row is for an amplification curve point.

plate Plate identifier. There is one identifier for each of the four data sets.

well Well identifier, i.e. the position within a PCR plate. This information was not available from the original publication, thus all values are NA.

target Target identifier. In all data sets the target is an amplicon consisting of soybean Lectin endogene "Le1".

dye Type of fluorescence dye, in this data set it is always SYBR Green I master mix (Roche) ("SYBR").

inhibitor Name of the molecule used as PCR inhibitor. In the case of the dilution series the value is "none".

inhibitor_conc Inhibitor concentration. Units are % (v/v) for isopropanol, and ug/mL for tannic acid.

sample Name of the biological sample. Samples have a simple consecutive identifier: S1, S2, ..., S5.

sample_type Sample type. All reactions are standard curves, i.e. "std".

replicate Replicate identifier.

copies Standard copy number of the amplicon.

dilution Dilution factor. Higher number means greater dilution, e.g. 5 means a 1:5 (five-fold) dilution relative to most concentrated standard.

cycle PCR cycle.

fluor Raw fluorescence values.

Source

[doi:10.1093/nar/gkr775](https://doi.org/10.1093/nar/gkr775)

Examples

lievens

Index

* **datasets**
 lievens, [2](#)

 lievens, [2](#)

 tibble, [3](#)